



## SARTORIUS

|            | eLine     |         | Picus, Picus NxT |               | ePet   |         | mLine                |               |
|------------|-----------|---------|------------------|---------------|--------|---------|----------------------|---------------|
| Channels ▶ | 1         | 8,12    | 1                | 4,8,12        | 1      | 8,12    | 1                    | 8,12          |
| Tip Type   | Tip Sizes |         |                  |               |        |         |                      |               |
| <b>A</b>   | 0.2-10    | 0.2-10  | 0.2-10           | 0.2-10        | 0.2-10 | 0.2-10  | 0.1-3, 0.5-10        | 0.5-10        |
| <b>U</b>   | 0.2-10    | 0.2-10  | 0.2-10           | 0.2-10        |        |         | 0.5-10               | 0.5-10        |
| <b>C</b>   | 5-120     | 5-120   | 5-120            | 5-120         | 5-100  | 5-100   | 2-20, 10-100, 20-200 | 5-100, 20-200 |
| <b>D</b>   |           |         |                  |               |        |         | 10-100, 20-200       |               |
| <b>E</b>   | 5-120     | 5-120   | 5-120            | 5-120         |        | 25-250  | 20-200               |               |
| <b>G</b>   | 10-300    | 5-120   | 10-300           | 5-120, 10-300 |        | 25-250  | 20-200               | 30-300        |
| <b>H</b>   | 5-120     | 5-120   | 5-120            | 5-120         | 5-100  |         | 10-100, 20-200       | 5-100         |
| <b>L</b>   | 50-1000   |         | 50-1000          |               |        |         | 100-1000             |               |
| <b>I</b>   | 50-1000   |         | 50-1000          |               |        |         |                      |               |
| <b>K</b>   |           | 50-1200 |                  | 50-1200       |        | 50-1200 |                      |               |
| <b>W</b>   |           |         |                  |               |        |         | 500-5000             |               |
| <b>X</b>   |           |         | 500-5000         |               |        |         | 500-5000             |               |
| <b>Y</b>   |           |         | 1000-10000       |               |        |         | 1000-10000           |               |

µL volume  
ranges shown  
**µL**



## SARTORIUS (Proline)

|            | Proline Electronic |        | Proline            |                | Proline Plus         |                |
|------------|--------------------|--------|--------------------|----------------|----------------------|----------------|
| Channels ▶ | 1                  | 8,12   | 1                  | 4,8,12         | 1                    | 8,12           |
| Tip Type   | Tip Sizes          |        |                    |                |                      |                |
| <b>A</b>   | 0.2-10             | 0.2-10 | 0.1-2.5, 0.5-10    | 0.5-10         | 0.1-3, 0.5-10        | 0.5-10         |
| <b>U</b>   |                    |        |                    |                | 0.1-3, 0.5-10        | 0.5-10         |
| <b>C</b>   | 5-100              | 5-100  | 2-20, 5-50, 10-100 | 5-50           | 2-20, 10-100, 20-200 | 10-100, 20-200 |
| <b>D</b>   |                    |        | 10-100, 20-200     |                | 10-100, 20-200       |                |
| <b>E</b>   | 5-100              | 5-100  | 20-200             |                | 20-200               |                |
| <b>F</b>   |                    |        | 20-200             |                |                      | 5-250          |
| <b>G</b>   | 25-250             | 10-250 | 50-200, 50-250     | 50-250, 50-300 | 20-200, 20-300       | 50-250, 30-300 |
| <b>H</b>   | 5-100              |        | 5-50, 50-200       |                | 2-20, 10-100         | 10-100         |
| <b>L</b>   | 50-1000            |        |                    |                | 100-1000             |                |
| <b>I</b>   | 50-1000            |        | 100-1000           |                | 100-1000             |                |
| <b>K</b>   | 50-1200            |        |                    |                |                      |                |
| <b>W</b>   | 100-5000           |        | 100-5000           |                |                      |                |
| <b>X</b>   |                    |        |                    |                | 500-5000             |                |
| <b>Y</b>   |                    |        |                    |                | 1000-10000           |                |



## BRAND

|            | Transferpette              |                                    | Electronic   |               | Transferpette S |                      |
|------------|----------------------------|------------------------------------|--------------|---------------|-----------------|----------------------|
| Channels ▶ | 1                          | 8,12                               | 1            | 4, 8,12       | 1               | 8, 12                |
| Tip Type   | Tip Sizes                  |                                    |              |               |                 |                      |
| <b>A</b>   | 0.1-1, 0.1-10              | 0.5-10                             |              | 0.5-10        | 0.1-1, 0.5-10   | 0.5-10               |
| <b>B</b>   | 0.5-10, 2-20               | 2-20                               | 0.5-10, 2-20 |               | 2-20            |                      |
| <b>C</b>   | 2-20, 5-50, 10-100, 20-200 | 2-20, 2.5-25, 5-50, 10-100, 20-200 | 20-200       | 5-100, 10-200 | 10-100, 20-200  | 5-50, 10-100, 20-200 |
| <b>D</b>   |                            |                                    | 20-200       | 5-100, 10-200 | 10-100, 20-200  |                      |
| <b>E</b>   | 20-200                     | 20-200                             | 20-200       |               | 20-200          | 20-200               |
| <b>F</b>   | 20-200, 25-250             |                                    |              |               |                 |                      |
| <b>G</b>   | 20-200, 25-250             | 30-300                             |              | 15-300        |                 |                      |
| <b>H</b>   | 5-50, 10-100, 20-200       | 5-50, 10-100, 20-200               | 20-200       | 5-100, 10-200 | 20-200          | 10-100, 20-200       |
| <b>L</b>   | 100-1000                   |                                    | 100-1000     |               | 100-1000        | 100-1000             |
| <b>I</b>   | 100-1000                   |                                    |              |               | 100-1000        | 100-1000             |



## THERMO /Finnpipette

|          | Digital, F1, F2, F3                |            |      | Focus                    |        | Novus              |        |      |
|----------|------------------------------------|------------|------|--------------------------|--------|--------------------|--------|------|
|          | 1                                  | 8,12       | 16   | 1                        | 8,12   | 1                  | 8,12   | 16   |
| Tip Type | Tip Sizes                          |            |      |                          |        |                    |        |      |
| <b>A</b> | 0.2-2, 0.5-10                      | 0.5-10*    |      | 0.3-3, 0.5-5, 1-10       | 1-10   | 1-10               | 1-10   |      |
| <b>C</b> | 0.5-10, 2-20, 5-50, 10-100, 20-200 | 5-50, 1-10 |      | 1-10, 3-30, 5-50, 10-100 | 5-50   | 1-10, 5-50, 10-100 | 5-50   |      |
| <b>E</b> | 2-20, 5-50, 10-100, 20-200         | 5-50       |      |                          |        |                    |        |      |
| <b>G</b> |                                    | 50-300     |      | 30-300                   | 30-300 | 30-300             | 30-300 |      |
| <b>H</b> | 2-20, 5-50, 10-100, 20-200         | 5-50       |      | 1-10, 3-30, 5-50, 10-100 |        | 10-100             | 5-50   |      |
| <b>L</b> |                                    |            |      | 100-1000                 |        | 100-1000           |        |      |
| <b>I</b> | 100-1000, 200-1000                 |            |      | 100-1000                 |        | 100-1000           |        |      |
| <b>T</b> |                                    |            | 1-50 |                          |        |                    |        | 1-50 |

\* for digital only



# EPENDORF

|            | Electronic    | Reference, Reference 2 |        | Research Pro     |        | Research, Research Plus, Xplorer, Xplorer Plus |                |
|------------|---------------|------------------------|--------|------------------|--------|--|----------------|
| Channels ▶ | 1             | 1                      | 8, 12  | 1                | 8, 12  | 1  | 8, 12          |
| Tip Type   | Tip Sizes     |                        |        |                  |        |  |                |
| <b>A</b>   | 0.5-10        | 0.1-2.5, 0.5-10        | 0.5-10 | 0.5-10           | 0.5-10 | 0.1-2.5, 0.5-10                                | 0.5-10         |
| <b>B</b>   | 0.5-20        | 2-20                   |        | 0.5-10           | 0.5-10 |  |                |
| <b>C</b>   | 5-100, 20-200 | 2-20, 10-100, 50-200   | 10-100 | 5-100            | 5-100  | 2-20, 5-100, 10-100, 20-200                    | 5-100, 10-100  |
| <b>D</b>   | 5-100, 20-200 | 2-20, 10-100, 50-200   | 10-100 | 5-100            | 5-100  | 2-20, 5-100, 10-100, 20-200                    | 10-100         |
| <b>E</b>   | 5-100, 20-200 | 2-20, 10-100, 50-200   | 10-100 | 5-100            | 5-100  | 2-20, 5-100, 10-100, 20-200                    | 10-100         |
| <b>F</b>   | 20-200        |                        |        |                  |        | 10-100, 20-200                                 |                |
| <b>G</b>   | 20-200        | 50-200, 30-300         | 30-300 | 20-300           | 20-300 |  | 15-300, 30-300 |
| <b>H</b>   | 5-100, 20-200 | 2-20, 10-100, 5-200    |        | 5-100            | 5-100  | 2-20, 10-100, 20-200                           | 10-100         |
| <b>L</b>   |               | 100-1000               |        | 50-1000          |        | 100-1000                                       |                |
| <b>I</b>   | 50-1000       | 100-1000               |        | 50-1000, 50-1200 |        | 100-1000                                       |                |
| <b>W</b>   |               |                        |        | 100-5000         |        | 500-5000                                       |                |

µL volume  
ranges shown  
**µL**



# GILSON

|            | Pipetman Classic, G, L           |              | Pipetman Ultra          |        | Neo                     | Concept |        | Pipetman M |              |
|------------|----------------------------------|--------------|-------------------------|--------|-------------------------|---------|--------|------------|--------------|
| Channels ▶ | 1                                | 8            | 1                       | 8, 12  | 1                       | 1       | 8, 12  | 1          | 8            |
| Tip Type   | Tip Sizes                        |              |                         |        |                         |         |        |            |              |
| <b>A</b>   | 0.2-2, 1-10                      |              | 0.2-2, 1-10             |        | 0.2-2, 1-10             | 0.5-10  | 1-10   | 0.5-10     |              |
| <b>U</b>   | 0.2-2, 1-10                      |              | 0.2-2, 1-10             |        | 0.2-2, 1-10             |         |        |            |              |
| <b>C</b>   | 2-20, 10-200,<br>20-200, 50-200  | 1-20, 20-200 | 2-20, 20-100,<br>20-200 | 1-20   | 2-20, 20-100,<br>20-200 | 5-100   | 10-100 | 2-20       | 1-20, 20-200 |
| <b>E</b>   | 50-200                           |              | 50-200                  |        |                         |         |        |            |              |
| <b>G</b>   | 50-200                           |              | 20-200                  | 20-300 |                         | 20-300  | 30-300 |            |              |
| <b>H</b>   | 10-100, 20-100,<br>5-200, 50-200 |              | 20-100,<br>20-200       |        | 20-200                  | 5-100   | 10-100 |            |              |
| <b>L</b>   | 20-1000                          |              | 20-1000                 |        | 100-1000                |         |        | 100-1000   |              |
| <b>I</b>   | 20-1000                          |              | 20-1000                 |        | 100-1000                |         |        | 100-1000   |              |
| <b>O</b>   | 1000-5000                        |              |                         |        |                         |         |        |            |              |
| <b>X</b>   | 1000-5000                        |              |                         |        |                         |         |        |            |              |
| <b>Y</b>   | 1000-10000                       |              | 1000-10000              |        |                         |         |        |            |              |



## METTLER TOLEDO RAININ

|            | E-Man Hybrid         | EDP3 Plus | EDP3 with LTS |               | E4XLS        |              | EDP1                 |
|------------|----------------------|-----------|---------------|---------------|--------------|--------------|----------------------|
| Channels ▶ | 1                    | 1         | 1             | 8, 12, 16, 24 | 8            | 12           | 1                    |
| Tip Type   | Tip Sizes            |           |               |               |              |              |                      |
| <b>A</b>   | 0.5-10               |           |               |               |              |              | 0.5-10               |
| <b>Q</b>   | 0.1-2, 0.5-10, 2-20  |           |               | 0.5-10, 2-20  | 0.5-10       | 0.5-10, 2-20 | 0.5-10               |
| <b>C</b>   | 2-20, 10-100, 20-200 |           |               |               |              |              | 2-20, 10-100, 20-200 |
| <b>E</b>   | 20-200               |           |               |               |              |              | 20-200               |
| <b>R</b>   | 10-100, 20-200       |           |               | 5-50, 20-200  | 5-50, 20-200 | 5-50, 20-200 | 5-50, 20-200         |
| <b>F</b>   | 20-200               |           |               |               |              |              |                      |
| <b>G</b>   | 20-200, 20-300       |           |               |               |              |              | 20-300               |
| <b>H</b>   | 2-20, 10-100, 20-200 |           |               |               |              |              | 2-20, 10-100, 20-200 |
| <b>L</b>   | 100-1000             |           |               |               |              |              | 100-1000             |
| <b>I</b>   | 100-1000             |           |               |               |              |              | 100-1000             |
| <b>S</b>   | 100-1000             |           |               | 100-1000      | 100-1000     | 100-1000     | 100-1000             |
| <b>V</b>   |                      |           |               |               | 20-300       | 20-300       |                      |

μL volume  
ranges shown  
**μL**



## METTLER TOLEDO RAININ

|            | EDP1 with LTS  | Pipet Lite SL | Pipet Lite/LTS       | XLS+     | Classic            |
|------------|----------------|---------------|----------------------|----------|--------------------|
| Channels ▶ | 1              | 8, 12, 18, 24 | 1                    | 1, 8, 12 | 1, 8, 12           |
| Tip Type   | Tip Sizes      |               |                      |          |                    |
| <b>A</b>   |                |               | 0.1-2, 0.5-10        |          |                    |
| <b>Q</b>   | 0.5-10, 2-20   | 1-10, 2-20    |                      |          | 0.1-2, 0.5-10      |
| <b>C</b>   |                |               | 2-20, 10-100, 20-200 |          |                    |
| <b>E</b>   |                |               | 20-200               |          |                    |
| <b>R</b>   | 10-100, 20-200 | 5-50, 20-200  |                      |          | 0.1-2, 0.5-10      |
| <b>F</b>   |                |               | 20-200               |          |                    |
| <b>H</b>   |                |               | 2-20, 10-100, 20-200 |          |                    |
| <b>L</b>   |                |               | 100-1000             |          |                    |
| <b>I</b>   |                |               | 100-1000             |          |                    |
| <b>O</b>   |                |               |                      |          | 1000-5000          |
| <b>S</b>   | 100-1000       |               |                      | 100-1000 | 100-1000, 100-1200 |
| <b>X</b>   |                |               |                      |          | 1000-5000          |
| <b>Y</b>   |                |               |                      |          | 5000-10000         |
| <b>V</b>   |                |               | 20-300               | 20-300   |                    |



## SOCOREX

|            | Acura 825, Manual XS             |              | Acura 835 | Calibra 822                      | Calibra 825    | Electro 925    | Electro 955 |
|------------|----------------------------------|--------------|-----------|----------------------------------|----------------|----------------|-------------|
| Channels ▶ | 1                                | 8, 12        | 1         | 1                                | 8, 12          | 1              | 8, 12       |
| Tip Type   | Tip Sizes                        |              |           |                                  |                |                |             |
| <b>A</b>   | 0.1-2, 0.5-10                    |              |           |                                  |                | 0.5-10         | 0.5-10      |
| <b>B</b>   |                                  | 0.5-10       |           | 0.2-2, 1-10                      | 1-10           |                |             |
| <b>C</b>   | 1-10, 2-20, 5-50, 10-100, 20-200 | 5-50, 20-200 |           | 1-10, 2-20, 5-50, 10-100, 20-200 | 10-100, 20-200 | 2.5-50, 10-100 | 2.5-50      |
| <b>E</b>   | 20-200                           |              |           | 20-200                           | 20-200         | 10-200         |             |
| <b>F</b>   | 20-200                           | 20-200       |           |                                  |                |                |             |
| <b>G</b>   |                                  | 40-350       |           | 20-200                           | 20-200         |                |             |
| <b>H</b>   | 20-200                           | 20-200       |           | 20-200                           | 20-200         |                | 2.5-50      |
| <b>L</b>   | 100-1000                         |              |           | 100-1000                         |                |                |             |
| <b>I</b>   | 100-1000                         |              |           | 100-1000                         |                |                |             |
| <b>W</b>   |                                  |              | 500-5000  |                                  |                |                |             |
| <b>X</b>   |                                  |              | 500-5000  |                                  |                |                |             |



## HTL and LABNET

|            | Excel     |          | Labpette Discovery         |              | Discovery Pro |
|------------|-----------|----------|----------------------------|--------------|---------------|
| Channels ▶ | 1         | 8, 12    | 1                          | 8, 12        | 1, 8, 12      |
| Tip Type   | Tip Sizes |          |                            |              |               |
| <b>A</b>   | 0.5-10    | 0.5-10   | 0.2-2, 0.5-10              | 1-10         | 0.2-2, 0.5-10 |
| <b>U</b>   | 2-20      | 2-20     |                            |              |               |
| <b>C</b>   | 10-200    | 10-200   | 2-20, 5-50, 10-100, 20-200 | 5-50, 20-200 | 2-20, 20-200  |
| <b>E</b>   | 10-200    | 10-200   |                            |              |               |
| <b>G</b>   |           |          | 50-250                     | 50-300       |               |
| <b>H</b>   | 10-200    | 10-200   | 10-100, 20-200             | 20-200       |               |
| <b>L</b>   |           |          | 100-1000                   |              | 100-1000      |
| <b>I</b>   | 100-1200  | 100-1200 | 100-1000                   |              | 100-1000      |
| <b>O</b>   |           |          | 1000-5000                  |              |               |
| <b>Y</b>   |           |          | 1000-10000                 |              |               |

µL volume  
ranges shown  
**µL**



## HAMILTON

|            | Softgrip       |        | Softouch Electronic |        |
|------------|----------------|--------|---------------------|--------|
| Channels ▶ | 1              | 8, 12  | 1                   | 8, 12  |
| Tip Type   | Tip Sizes      |        |                     |        |
| <b>A</b>   | 0.2-2, 1-10    |        | 0.2-10              | 0.2-10 |
| <b>C</b>   | 2.5-25, 10-100 | 5-50   | 5-100               |        |
| <b>E</b>   | 10-100         |        | 5-100               |        |
| <b>F</b>   | 10-100         |        |                     |        |
| <b>G</b>   |                | 30-300 |                     | 25-250 |
| <b>H</b>   | 2.5-25, 10-100 | 5-50   | 5-100               |        |
| <b>L</b>   | 100-1000       |        | 50-1000             |        |
| <b>I</b>   | 100-1000       |        | 50-1000             |        |



## NICHIRYO

|                 | EX, EX Plus          | Le                   | Nichipet 7000  |
|-----------------|----------------------|----------------------|----------------|
| Channels ▶      | 1                    | 1                    | 8, 12          |
| <b>Tip Type</b> | <b>Tip Sizes</b>     |                      |                |
| <b>A</b>        | 0.2-2, 0.5-10        | 0.5-10               | 0.5-10         |
| <b>B</b>        | 0.5-10               |                      |                |
| <b>C</b>        | 2-20, 10-100, 20-200 | 2-20, 10-100, 20-200 | 5-50, 40-200   |
| <b>E</b>        | 20-200               | 20-200               |                |
| <b>F</b>        | 20-200               | 20-200               |                |
| <b>G</b>        | 20-200               | 20-200               | 40-200, 50-300 |
| <b>H</b>        | 10-100, 20-200       | 10-100, 20-200       | 5-50, 40-200   |
| <b>L</b>        | 100-1000             | 100-1000             |                |
| <b>I</b>        | 100-1000             | 100-1000             |                |
| <b>X</b>        | 1000-5000            |                      |                |



## VISTALAB

|                 | MLA              | Ovation             | Ovation Electronic |              |
|-----------------|------------------|---------------------|--------------------|--------------|
| Channels ▶      | 1                | 1                   | 1                  | 8, 12        |
| <b>Tip Type</b> | <b>Tip Sizes</b> |                     |                    |              |
| <b>A</b>        |                  | 1-10                | 0.5-20             | 0.5-20       |
| <b>U</b>        |                  | 0.2-2, 1-10         | 0.2-2              | 0.2-2        |
| <b>C</b>        | 10-50, 50-200    | 10-100, 20-200      | 2-125, 5-250       | 2-125, 5-250 |
| <b>E</b>        | 10-50, 50-200    | 20-200              | 2-125, 5-250       | 2-125, 5-250 |
| <b>F</b>        | 50-200           |                     |                    |              |
| <b>G</b>        |                  | 20-200              | 5-250              | 5-250        |
| <b>H</b>        | 10-50, 50-200    | 10-100, 20-200      | 2-125              | 2-125        |
| <b>L</b>        |                  | 100-1000, 201-1000* |                    |              |
| <b>I</b>        |                  | 100-1000*           | 25-1250*           |              |
| <b>O</b>        | 1000-5000        |                     | 100-5000           |              |
| <b>W</b>        | 1000-5000        |                     | 100-5000           |              |

\* requires Type A Nozzle Vistalab #9057-1016 100-1000µL, 201-1000µL

µL volume  
ranges shown  
**µL**



## DRAGONLAB

TopPette, MicroPette, MicroPette Plus

| Channels ▶      | 1, 8, 12                           |
|-----------------|------------------------------------|
| <b>Tip Type</b> | <b>Tip Sizes</b>                   |
| <b>A</b>        | 0.1-2.5, 0.5-10                    |
| <b>C</b>        | 2-20, 5-50, 10-100, 20-200, 50-200 |
| <b>E</b>        | 20-200, 50-200                     |
| <b>F</b>        | 20-200                             |
| <b>H</b>        | 20-200                             |
| <b>X</b>        | 1000-5000*                         |

\* Single channel only



## ACCUMAX

Smart, Accumax Pro, A, VA-FA

| Channels ▶      | 1, 8, 12                           |
|-----------------|------------------------------------|
| <b>Tip Type</b> | <b>Tip Sizes</b>                   |
| <b>A</b>        | 0.1-2.5, 0.5-10                    |
| <b>C</b>        | 2-20, 5-50, 10-100, 20-200, 50-200 |
| <b>E</b>        | 20-200, 50-200                     |
| <b>F</b>        | 20-200                             |
| <b>H</b>        | 20-200                             |



## POULTEN

Volac Ultra Variable

| Channels ▶ | 1                    |
|------------|----------------------|
| Tip Type   | Tip Sizes            |
| <b>B</b>   | 1-10                 |
| <b>C</b>   | 5-50, 20-100, 20-200 |
| <b>E</b>   | 20-200               |
| <b>F</b>   | 20-200               |
| <b>H</b>   | 20-200               |
| <b>I</b>   | 100-1000             |



## CAPP

Aero, Aero Comfort, Bravo, Eco

| Channels ▶ | 1                          | 8, 12         |
|------------|----------------------------|---------------|
| Tip Type   | Tip Sizes                  |               |
| <b>A</b>   | 0.1-2, 0.5-10              | 0.1-2, 0.5-10 |
| <b>B</b>   | 0.5-10                     |               |
| <b>C</b>   | 2-20, 5-50, 10-100, 20-200 | 5-50, 20-200  |
| <b>E</b>   | 20-200                     | 20-200        |
| <b>G</b>   | 50-300                     | 30-300        |
| <b>H</b>   | 20-200                     | 20-200        |
| <b>L</b>   | 100-1000                   |               |
| <b>I</b>   | 100-1000                   |               |



## KENDALL

Oxford Benchmate

| Channels ▶ | 1                      |
|------------|------------------------|
| Tip Type   | Tip Sizes              |
| <b>A</b>   | 0.1-2, 0.5-10          |
| <b>B</b>   | 0.5-10                 |
| <b>C</b>   | 10-100, 20-200, 40-200 |
| <b>E</b>   | 20-200                 |
| <b>F</b>   | 20-200                 |
| <b>G</b>   | 20-200, 40-200         |
| <b>H</b>   | 10-100, 20-200         |
| <b>I</b>   | 200-1000               |
| <b>O</b>   | 1000-5000              |



## CORNING

8/12-Pette Octapette Lambda

| Channels ▶ | 8, 12     | 8             | 1                    |
|------------|-----------|---------------|----------------------|
| Tip Type   | Tip Sizes |               |                      |
| <b>A</b>   |           |               | 0.1-2.5, 0.5-10      |
| <b>C</b>   | 20-200    | 25,50,100,200 | 2-20, 10-100, 20-200 |
| <b>E</b>   | 20-200    | 25,50,100,200 | 2-20, 10-100, 20-200 |
| <b>F</b>   | 20-200    |               | 10-100, 20-200       |
| <b>G</b>   | 20-200    |               |                      |
| <b>H</b>   | 20-200    |               | 2-20, 10-100, 20-200 |
| <b>I</b>   |           |               | 100-1000             |



## METTLER

Multimate

| Channels ▶ | 1                    |
|------------|----------------------|
| Tip Type   | Tip Sizes            |
| <b>A</b>   | 0.1-2.5, 0.5-10      |
| <b>C</b>   | 2-20, 10-100, 20-200 |
| <b>D</b>   | 2-20, 10-100, 20-200 |
| <b>E</b>   | 2-20, 10-100, 20-200 |
| <b>F</b>   | 10-100, 20-200       |
| <b>H</b>   | 2-20, 10-100, 20-200 |
| <b>I</b>   | 100-1000             |



## THERMO/Matrix

Impact 2 Shorty

| Channels ▶ | 1, 8         | 1            |
|------------|--------------|--------------|
| Tip Type   | Tip Sizes    |              |
| <b>U</b>   | 12.5, 30**   | 12.5, 30**   |
| <b>C</b>   | 2-125, 5-250 | 2-125, 5-250 |
| <b>E</b>   | 5-250        | 5-250        |
| <b>F</b>   | 5-250        | 5-250        |
| <b>G</b>   | 5-250        | 5-250        |
| <b>H</b>   | 2-125        | 2-125        |
| <b>J</b>   | 15-1250      | 15-1250      |

\*\* Not for 16 Channel Equalizer