


Robotic Pipet Tips for Workstations

- For Beckman® Workstations
- Automation Certified
- RNase/DNase Free
- Non pyrogenic
- Pre-tested medical grade resin
- Color coded racks





Robotic Pipet Tips

Labcon automation pipet tips (APT™) for Beckman® FX and NX are designed and tested to fit your workstation. Made from pure pre-tested medical grade resins to ensure the integrity of your work. Each lot is tested and certified free of RNase/DNase (Nuclease), and Endotoxins (non pyrogenic). Resins are pre-tested with a U.S.P. test method to ensure they are metal free, and they are free of Oleamide and Phalates. Unlike tips sold as “pre-sterile” these tips use a validated medical grade sterilization process. So you can use these tips in diagnostics, regulated testing, and medical screening procedures.

Labcon Robotic Pipet Tips are made of virgin polypropylene and furnished pre-racked in extra heavy automation grade polypropylene racks. Autoclavable.

- For Beckman® FX and NX
- Non pyrogenic, Endotoxin, Nuclease, Protease, Metal Free
- Pure medical grade resins
- Sterile (validated) and Non-sterile available


| Code | Item | Package | Qty/Pk | Sterile | Qty/Cs |
|------|----------|--|---------------------------------------|---------|--------|
| A | 8005-690 | P30 AP384 Style, 30µL capacity | 384 Tips/Rack, 10 Racks/Pkg, 5 Pkg/Cs | | 19,200 |
| A | 8005-695 | P30 AP384 Style, 30µL capacity | 384 Tips/Rack, 10 Racks/Pkg, 5 Pkg/Cs | ✓ | 19,200 |
| B | 8014-690 | P20 AP96 Style, 100µL capacity | 96 Tips/Rack, 10 Racks/Pkg, 5 Pkg/Cs | | 4,800 |
| B | 8014-695 | P20 AP96 Style, 100µL capacity | 96 Tips/Rack, 10 Racks/Pkg, 5 Pkg/Cs | ✓ | 4,800 |
| B | 8114-695 | P20 AP96 Style, with Filter, 20µL capacity | 96 Tips/Rack, 10 Racks/Pkg, 5 Pkg/Cs | ✓ | 4,800 |
| C | 8024-690 | P50XL AP96 Style, 90µL capacity | 96 Tips/Rack, 10 Racks/Pkg, 5 Pkg/Cs | | 4,800 |
| C | 8024-695 | P50XL AP96 Style, 90µL capacity | 96 Tips/Rack, 10 Racks/Pkg, 5 Pkg/Cs | ✓ | 4,800 |
| D | 8017-690 | P250XL AP96 Style, 220µL capacity | 96 Tips/Rack, 10 Racks/Pkg, 5 Pkg/Cs | | 4,800 |
| D | 8017-695 | P250XL AP96 Style, 220µL capacity | 96 Tips/Rack, 10 Racks/Pkg, 5 Pkg/Cs | ✓ | 4,800 |
| D | 8117-695 | P250XL AP96 Style, with Filter, 170µL capacity | 96 Tips/Rack, 10 Racks/Pkg, 5 Pkg/Cs | ✓ | 4,800 |



Every lot of APT™ tips are meticulously tested and inspected.

TIPS ARE SHOWN ACTUAL SIZE.
SEE PAGE 6 FOR DEFINITIONS OF ALPHABETICAL SYMBOLS THAT APPEAR BELOW THIS BOX.

Packaging Options



96 Rack