

culture tubes

SuperClear™ culture tubes and caps



These precision-molded, conical and round bottom tubes are perfect for a multitude of general laboratory uses. Precision crafted of extra clear virgin polypropylene in natural and colors or clear virgin polystyrene. Tube caps are made of low density polyethylene for a secure seal and are available in a fluted plug style or dual position design that permits sample aspiration.

Polypropylene tubes are autoclavable to 122°C and freezable to -80°C. Polystyrene tubes are not autoclavable but have a nominal temperature range between 90°C and -40°C. Polyethylene caps are not autoclavable, but are available pre-sterilized with many of the tube packages.

Polypropylene tubes can be centrifuged to 3000 RCF and polystyrene tubes will withstand 1500 RCF. (RCF; Relative Centrifugal Force, is not equivalent to RPM).

Non-pyrogenic



12x75mm Round Bottom and Conical tubes with dual position or plug style caps.



12 x 75mm, Conical

12 x 75mm, Round

Cap Options:



Plug Type



Dual Position

12x75mm Conical Tubes, 4mL Capacity, without Caps

Item	Description	Qty/Pk	Sterile	Qty/Cs
3360-800	Polystyrene, Clear, 1000/Bag	1000		1000
3361-800	Polypropylene, Natural, 1000/Bag	1000		1000
3385-335	Polystyrene, Clear, 25/Bag	25	✓	500
3386-335	Polypropylene, Natural, 25/Bag	25	✓	500

12x75mm Round Bottom Tubes, 5mL Capacity, without Caps

Item	Description	Qty/Pk	Sterile	Qty/Cs
3325-350	Polystyrene, Clear, 125/Bag	125		1000
3325-355	Polystyrene, Natural, 125/Bag	125	✓	1000
3326-350	Polypropylene, Clear, 125/Bag	125		1000
3326-355	Polypropylene, Natural, 125/Bag	125	✓	1000
3305-800	Polystyrene, Clear, 3000/Bag	3000		3000
3306-800	Polypropylene, Natural, 3000/Bag	3000		3000
3310-800	Polystyrene, Clear, 1000/Bag	1000		1000
3311-800	Polypropylene, Natural, 1000/Bag	1000		1000
3311-802	Polypropylene, Yellow, 1000/Bag	1000		1000
3311-803	Polypropylene, Orange, 1000/Bag	1000		1000
3311-804	Polypropylene, Red, 1000/Bag	1000		1000
3311-806	Polypropylene, Purple, 1000/Bag	1000		1000
3311-807	Polypropylene, Blue, 1000/Bag	1000		1000
3311-808	Polypropylene, Green, 1000/Bag	1000		1000

12x75mm Round Bottom Tubes, 5mL Capacity, with Separate Dual-Position Caps, Sterile

Item	Description	Qty/Pk	Sterile	Qty/Cs
3315-800	Polystyrene, Clear, 1000/Bag	1000		1000
3316-800	Polypropylene, Natural, 1000/Bag	1000		1000

12x75mm Round Bottom Tubes, 5mL Capacity, attached Dual-Position Caps, Sterile

Item	Description	Qty/Pk	Sterile	Qty/Cs
3328-355	Polystyrene, Clear, 125/Bag	125	✓	1000
3329-355	Polypropylene, Natural, 125/Bag	125	✓	1000
3335-335	Polystyrene, Clear, 25/Bag	25	✓	500
3336-335	Polypropylene, Natural, 25/Bag	25	✓	500

12x75mm Round Bottom Tubes, 5mL Capacity, attached Plug-Type Caps, Sterile


Item	Description	Qty/Pk	Sterile	Qty/Cs
3342-355	Polystyrene, Clear, 25/Bag	25	✓	500
3343-355	Polypropylene, Natural, 25/Bag	25	✓	500

Dual-Position Caps for 12x75mm Tubes

Item	Description	Qty/Pk	Sterile	Qty/Cs
3300-800	Polyethylene, Clear, 1000/Bag	1000		3000



Labcon SuperClear™ tubes are clear enough to see color gradients. Polypropylene tubes available in colors.



13x100mm and 17x100mm Round Bottom tubes with dual position and plug-type caps.



Heavy walls on Labcon tubes protect your sample



Other brands of culture tubes have thin walls that crush easily

- Ultra clear resin
- Colored tubes available
- Medical grade (U.S.P.) packaging
- Thicker walls will not crush

13x100mm Tubes, 8mL Capacity, without Caps

Item	Description	Qty/Pk	Sterile	Qty/Cs
3410-800	Polystyrene, Clear, 1000/Bag	1000		1000
3411-800	Polypropylene, Natural, 1000/Bag	1000		1000
3425-350	Polystyrene, Clear, 125/Bag	125		1000
3426-350	Polypropylene, Natural, 125/Bag	125		1000
3425-355	Polystyrene, Clear, 125/Bag	125	✓	1000
3426-355	Polypropylene, Natural, 125/Bag	125	✓	1000

17x100mm Tubes, 14mL Capacity, without Caps

Item	Description	Qty/Pk	Sterile	Qty/Cs
3510-800	Polystyrene, Clear, 1000/Bag	1000		1000
3511-800	Polypropylene, Natural, 1000/Bag	1000		1000
3525-350	Polystyrene, Clear, 125/Bag	125	✓	1000
3526-350	Polypropylene, Natural, 125/Bag	125	✓	1000

17x100mm Tubes, 14mL Capacity, with Polyethylene Dual-Position Caps

Item	Description	Qty/Pk	Sterile	Qty/Cs
3515-800	Polystyrene, Clear, 1000/Bag	1000		1000
3516-800	Polypropylene, Natural, 1000/Bag	1000		1000
3528-355	Polystyrene, Clear, 125/Bag	125	✓	1000
3529-355	Polypropylene, Natural, 125/Bag	125	✓	1000
3535-335	Polystyrene, Clear, 25/Bag	25	✓	500
3536-335	Polypropylene, Natural, 25/Bag	25	✓	500

Plug-Type Caps for 13x100mm and 12x75mm Tubes, Polyethylene


Item	Description	Qty/Pk	Sterile	Qty/Cs
3302-800	Natural, 1000/Bag	1000		3000
3302-801	White, 1000/Bag	1000		3000
3302-802	Yellow, 1000/Bag	1000		3000
3302-803	Orange, 100/Bag	1000		3000
3302-804	Red, 1000/Bag	1000		3000
3302-806	Blue, 1000/Bag	1000		3000
3302-807	Purple, 1000/Bag	1000		3000
3302-808	Green, 1000/Bag	1000		3000
3302-870	Assorted, 500/Bag	500		3000

Dual-Position Caps for 17x100mm Tubes

Item	Description	Qty/Pk	Qty/Cs
3500-800	Polyethylene, Natural, 1000/Bag	1000	1000

Plug-Type Caps for 17x100mm Tubes

Item	Description	Qty/Pk	Qty/Cs
3055-870	Polyethylene, Natural, 500/Bag	500	1000



17x100mm Round Bottom tube with dual position caps.



Colored plug-type caps for 13x100mm and 12x75mm tubes