



<1PPB

Metal Free* Centrifuge Tubes with ViewPoint™ labeling surface

Our Metal Free tubes are made with ultra clean resins for the tube and cap and processed to ensure that common trace metals are kept below an ICP-MS detection level of one part per billion. These tubes are perfect for storage of samples and for applications requiring the assurance of low background metals. Available in sterile racks or bags. Metals included in testing are Calcium, Magnesium, Zinc, Iron, Manganese, Copper, Aluminum, Silicon, Nickel, Vanadium, Sodium, Phosphorus, Cobalt, Chromium, Potassium, Lithium, and Lead. Autoclavable (122°C) and freezable (-80°C).

Spins to: 20,000 RCF

Item	Package	Qty/Pk	Sterile	Qty/Cs
50mL Metal Free* Centrifuge Tubes with Flat Caps	3184-345 50/Bag, 10 Bags/Case	50	✓	500
	3194-335 25/Rack, 20 Racks/Case	50	✓	500

* Certified to less than 1PPB for major metals

Spins to: 17,000 RCF

Item	Package	Qty/Pk	Sterile	Qty/Cs
15mL Metal Free* Centrifuge Tubes with Flat Caps	3134-345 50/Bag, 10 Bags/Case	50	✓	500
	3144-335 25/Rack, 20 Racks/Case	50	✓	500

* Certified to less than 1PPB for major metals

