

centrifugation

centrifuge tubes for special applications

USP CFR END RD ATP 122° -80°



UVsafe™ Amber 15mL and 50mL Centrifuge tubes for Light Sensitive samples.

These tubes are designed to block 100% of U.V. light. They can be used instead of foil to shield samples sensitive to light radiation. They have been lab tested and certified to block U.V. light, but will still pass visible light so you can see your sample volume. Autoclavable (122°C) and Freezable (-80°C).



Spins to:
12,500
RCF

	Item	Package	Qty/Pk	Sterile	Qty/Cs
50mL UVsafe™ Amber Centrifuge Tubes with Plug Style Caps	3182-445	50/Bag, 10 Bags/Case	50	✓	500
	3192-435	25/Rack, 20 Racks/Case	50	✓	500
	3192-871	Bulk, Tubes & Caps in Separate Bags	500		500

	Item	Package	Qty/Pk	Sterile	Qty/Cs
15mL UVsafe™ Amber Centrifuge Tubes with Plug Style Caps	3132-445	50/Bag, 10 Bags/Case	50	✓	500
	3142-435	25/Rack, 20 Racks/Case	50	✓	500
	3142-871	Bulk, Tubes & Caps in Separate Bags	500		500